ACRYPOL HI

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
Fabbricante

Pericolo

Ruthinum® - Dental Manufacturing S.p.A.

Via Cà Mignola Nuova 1699 - 45021 Badia Polesine (RO) ITALY

Tel. +39 0425 51628 - Fax +39 0425 590156 - www.ruthinum.it - info@ruthinum.it

(GB) HEAT CURING ACRYLIC RESIN "HIGH IMPACT" WITH A HIGH MOLECULAR WEIGHTFOR COMPLETE, PARTIAL END COMBINED REMOVABLE DENTAL PROSTHESES AND RELINING.

DIRECTIONS: Heat curing acrylic resin liquid and powder "HIGH IMPACT" with a high molecular weight for total and/or partial removable prostheses with reduced thickness, and for combined prostheses and relinings.

CONTRAINDICATIONS:

- Direct contact of non-polymerised material in the intraoral area.
- In case of proven allergy to components of Acry Pol HI.
- All uses not indicated.

PREPARATION: Mixing ratio: 21g powder (30ml in volume) and 10ml liquid. Pour the desired quantity of liquid into a glass container, add the powder by tilting the container and tapping it on the support surface until a thin dry layer of powder remains on the surface. Wait for 20-30 seconds, remove any excess powder, mix for 30 seconds, cover and leave to rest. At a temperature of 23°C, the mixture has a rest time of 15-20 minutes, and remains plastic (processing time) for 10-15 minutes.

N.B.: the temperature has a considerable influence on the rest times and plasticity of the dough.

FLASK PREPARATION: After eliminating any traces of wax and washing thoroughly with boiling water, let the flask cool to a temperature between 20°C and a maximum of 40°C, and insulate the plaster with Acry Film or similar insulating material. Next, pack the prosthesis mould with an adequate quantity of mixture. The two parts of the flask are then joined together with one or two pressing operations (at least 2 minutes), using a hydraulic press at a pressure of 120/150 kg/cm².

RECOMMENDED POLYMERISATION: Place the flask in water at room temperature until completely immersed. Slowly bring the water to the boil in at least 45 minutes. Simmer for 30 minutes, then leave to cool slowly in the water for another 30 minutes. Then, removing it from the water, cool the muffle at room temperature.

FINISHING AND POLISHING: Finish as usual, with burs, rubber points, brushes and with pumice. Mirror polish with the help of brushes and polishing products (Acry Lux Crema, Acry lux Liquido, Acry Lux Pasta).

N.B.: Observing the above procedures will result in a denture with high stability and excellent tensile strength (on average, twice as much as a normal resin for removable prostheses).

Powder: Store in a cool and dry place

Liquid: Store at a maximum temperature of 30°C

The product is designed for professional use in the dental sector. Follow the instructions for use. The manufacturer shall accept no responsibility for damages arising from different or inadequate use. The user must personally check the suitability of the product for the intended use, particularly if these uses are not reported in the instructions for use. Do not use after expiration. Warning: Any serious accident that may occur with the device must be reported by the user to the manufacturer and the competent authority of the Member State where the user and/or patient are established.

LIQUID: Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing vapours. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. Store in a well-ventilated place. Keep cool. Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains methyl methacrylate, ethylene dimethacrylate.

POWDER: Do not breathe dust. IF ON SKIN: Wash with plenty of soap and water.

Contains benzoyl peroxide, methyl methacrylate. May produce an allergic reaction.

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Ruthinium® - Dental Manufacturing S.p.A.