Stargel HP Hard

A

+

B

SILICONE PER ADDIZIONE A+B 23Sh.A
A-SILICON A+B 23Sh.A

SILICONE PAR ADDITION A+B 23Sh.A

SILICONA DE ADICIÓN A+B 23Sh.A

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A+B-ADDITIONSSILIKON 23Sh.A

(I) Istruzioni d'uso

(GB) Instructions for Use

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(D) Gebrauchsanweisung

Ruthinium® - Dental Manufacturing S.p.A.

Via Cà Mignola Nuova 1699 - 45021 Badia Polesine (RO) ITALY

Tel. +39 0425 51628 - Fax +39 0425 590156 - www.ruthinium.it - info@ruthinium.it

Stargel HP Hard

A
+
B

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 Fabbricante

(GB) A-SILICON A+B 23Sh.A

PREPARATION: Place the model to be duplicated in the upper part of an open muffle.

DOSAGE: Pour equal quantities of white liquid "A" (Activator) and coloured liquid "B"(Base) into a container.

MIXING: Mix the two components slowly for 60 seconds until a uniform colour is obtained. Mixing too energetically may lead to the inclusion of air bubbles, altering the quality of the duplication.

FILLING: Begin to pour the mixture onto the model placed in a muffle through one of the holes; continue to pour the mixture, raising the cup in a perpendicular direction in respect of the muffle, up to a height of about 35 cm, in order to eliminate any possible microscopic bubbles.

DUPLICATE-MODEL SEPARATION: Wait until the perfect polymerisation of the product - about 30 minutes - before removing the model with a jet of compressed air directed between the model and the STARGEL HP HARD. Never use mechanical force with a spatula or pliers since this operation could alter the duplicate.

CASTING THE DUPLICATE: The casting of the duplicate is carried out with the impression placed in its muffle.

STARGEL HP HARD after polymerisation are compatible with all plasters, coatings or cast resins.

STORAGE AND CONSERVATION: STARGEL HP HARD can be conserved for 24 months from production date at a temperature between 5°C and 27°C.

IMPORTANT NOTES: Work on perfectly clean models, steam cleaned if possible, especially when oils have been used for carrying out milling operations. In some cases, oils may inhibit reaction with the catalyst. It is recommended to wait until the model to be duplicated is perfectly dry. If boxing has to be carried out, use boxing wax or paste: do not use silicone for condensation, in order to avoid any phenomenon of incompatibility between the two silicones (Adding / Condensation). In the case of the presence (on the model to be duplicated) or the use (in boxing) of a resin or of an auto or photo-polymerising composite, wait in any case until the perfect stabilisation of these products before placing STARGEL HP HARD in contact. Some of these products could inhibit the polymerising reaction and alter the quality of the material.

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