



WHAT DOES IT DO?

THE INNOVATIVE SYSTEM THAT ALLOWS YOU TO REPOSITION ANY 3D-PRINTED MODEL INTO YOUR TRADITIONAL ARTICULATOR USING A MAGNET-BASED SYSTEM.

OPTIMIZED FOR 3D PRINTING OF MODELS, IT OFFERS A **UNIQUE** SOLUTION ON THE MARKET.

THE SYSTEM CONSISTS OF TWO PARTS: **SOFTWARE** AND **HARDWARE**.

THE SOFTWARE PART CAN BE DOWNLOADED FROM: WWW.BLENDERFORDENTAL.COM



SCAN THE **QR CODE**AND DOWNLOAD THE **SOFTWARE**



IN THE **DACOS SOFTWARE**, USERS WILL FIND A BUTTON THAT WILL DIRECT THEM TO THE **RUTHINIUM E-COMMERCE** SITE FOR PURCHASING THE KIT.



HOW DOES IT WORK?

THE USE OF **DACOS** IS BASED ON THE INNOVATIVE USE OF THE **DACOS MODULE** AND THE **DACOS PLATE.**

THIS SYSTEM ALLOWS YOU TO GENERATE A 3D-PRINTED MODEL WITH SPECIFIC ATTACHMENTS IN WHICH MAGNETS WILL BE LOCATED, ALLOWING THE MODEL TO EASILY AND SECURELY ATTACH TO THE ARTICULATOR.

THE **DACOS PLATE** IS DESIGNED TO **PERFECTLY** REFLECT THE GEOMETRY AND COORDINATES SPECIFIED IN THE **DACOS MODULE**.

THIS TECHNOLOGY ENSURES **UNPRECEDENTED** PRECISION IN THE REPOSITIONING OF MODELS, FAITHFULLY REFLECTING THE **ORIGINAL DIGITAL DESIGN** WHILE MAINTAINING THE **SPATIAL ACCURACY** REQUIRED FOR HIGH-QUALITY INTERVENTIONS AND CHECKS.



FOR ALL ARTICULATORS

DACOS IS COMPATIBLE WITH ALL TYPES OF ARTICULATORS THANKS TO ITS SELF-INSTALLATION PROCESS, WHICH INVOLVES PLASTERING THE DACOS PLATE IN THE DESIRED ARTICULATOR. REQUIRED ONLY FOR THE FIRST USE.

DACOS OFFERS SEVERAL HEIGHT OPTIONS:

- DACOS 40MM: COMPATIBLE WITH ARTICULATORS THAT HAVE A MINIMUM Z HEIGHT OF 65MM.
- DACOS 60MM: COMPATIBLE WITH ARTICULATORS THAT HAVE A MINIMUM Z HEIGHT OF 90MM.

THIS ALLOWS USERS TO UTILIZE **DACOS** ON A WIDE RANGE OF ARTICULATORS, ENSURING **PRECISION** AND **EASE OF USE** WITH JUST ONE-TIME INSTALLATION.

INCLUDED ANATOMICAL PLANES

THANKS TO THE **CAMPER** AND **FRANKFURT** PLANES AND THE **PILOT CALIBRATOR, DACOS** CAN BE INSTALLED ACCORDING TO THE DESIRED OCCLUSAL PLANE.







PRINT YOUR MASTER CALIBRATOR

CUSTOMISED ACCURACY: ELIMINATES THE VARIABLES OF EXTERNAL PRINTING PROCESSES BY CREATING A CALIBRATOR PERFECTLY SUITED TO YOUR SYSTEM ENSURING THAT THE INSTALLATION IS ALWAYS **PRECISE AND CORRECT.**


FREE FILE: DOWNLOAD THE STL FILE AND PRINT IT DIRECTLY TO YOUR 3D PRINTER TO BEST SYNCHRONISE THE DACOS SYSTEM WITH YOUR CONFIGURATION.

FLEXIBLE OFFSETS: CHOOSE BETWEEN TIGHT, STANDARD OR LOOSE TO CALIBRATE TO YOUR NEEDS (STANDARD RECOMMENDED).

DEVELOPED FOR 3D PRINTING

DACOS HAS BEEN DEVELOPED WITH THE PRINCIPLE OF 'INTERFERENCE-FREE SUPPORT'.

METICULOUSLY DESIGNED TO INTEGRATE PERFECTLY WITH 3D PRINTING NESTING, INCLUDING THOSE THAT USE SUPPORTS, THE INTERIOR OF THE ATTACHMENTS IS FREE OF SUPPORTS IN CRITICAL AREAS AND FEATURES REDUCED HORIZONTAL EXTENSIONS BELOW THE DEFORMATION THRESHOLD.



UNIQUENESS

INNOVATION THAT SETS STANDARDS.

DACOS IS THE ONLY SYSTEM THAT FACILITATES THE OCCLUSAL CHECK OF DIGITALLY DESIGNED AND PRODUCED DENTURES, IN YOUR ARTICULATOR, WITH PRECISION AND EASE.

UNIVERSAL

DACOS IS NOT ONLY COMPATIBLE WITH ALL ARTICULATORS, BUT IT IS ALSO DESIGNED TO BE USED BY EVERY DENTAL CAD DESIGNER, REGARDLESS OF THE SOFTWARE. THE DEDICATED DACOS MODULE ON BLENDER FOR DENTAL ALLOWS YOU TO IMPORT YOUR MODELS, POSITION THEM WITH PRECISION, AND FINALIZE THE PROCESS WITH JUST A FEW CLICKS.

AFTER CLICKING 'FINISH UPPER' AND 'FINISH LOWER', THE MODELS ARE READY TO BE RE-EXPORTED INTO YOUR USUAL SOFTWARE.

THIS FLEXIBILITY ENSURES A SEAMLESS AND INTEGRATED WORKFLOW, **OPTIMIZING BOTH TIME AND RESOURCES.**



GLUE JET

GLUE JET IS THE TOOL THAT WILL GUIDE YOU STEP BY STEP DURING THE **MAGNET PLACEMENT PROCESS,** ENSURING THEY ARE ALWAYS ORIENTED CORRECTLY.

THANKS TO THIS DEVICE, THE USER CAN INSTALL THE **MAGNETS** IN THE CORRECT DIRECTION **WITHOUT WORRYING** ABOUT THE **REPULSION BETWEEN POLES**.

GLUEJET IS COMPOSED OF TWO MAIN ELEMENTS:

- GLUEJET DISPENSER INSERTS THE MAGNET INTO THE DESIGNATED POSITION WITHIN THE MODE!
- GLUEJET PLUNGER REFLECTS EXACTLY THE ORIENTATION OF THE MAGNET ALREADY PRESENT IN THE DACOS PLATE.

WHEN USING **GLUEJET**, THE USER EXTRACTS THE MAGNET THAT WILL BE PLACED, WHICH IS ATTRACTED BY THE MAGNET IN THE **PLUNGER**, **ENSURING THAT THE ORIENTATION IS ALWAYS CORRECT.**





info@ruthinium.it | +39 0425 51628 ruthinium | shop.ruthinium.com













