# Ruthinium **FIBRA**

## Ruthinium **FIBRA** **KIT**

---

Ruthinium  
**FIBRA****MAT**

Ruthinium  
**FIBRA****MAT**  
KIT

Ruthinium  
**FIBRA****UHM**

Ruthinium  
**FIBRA****UHM**  
KIT

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## (GB) INFORMATION

### **RUTHINIUM FIBRA UHM**

Screw-retained substructures for implant-supported Toronto Bridge and Natural Bridge.

### **RUTHINIUM FIBRA MAT**

Reinforcements for tissue-supported dentures.

#### CONTRAINDICATIONS

- Not indicated in case of proven allergy to the components of Ruthinium Fibra.
- All uses that are not indicated.

#### Description

Contents of Ruthinium Fibra Kit:

- 2 sheets of Ruthinium Fibra UHM (25cm x 20cm)
- 3 sheets of Ruthinium Fibra MAT (10cm x 10cm)
- 2 bottles of Ruthinium Fibra Base Resin (25ml)
- 1 bottle of Ruthinium Fibra Catalyst (25ml)
- 1 silicone spatula
- 1 roller

#### Package contents:

- Ruthinium Fibra UHM - 3 sheets (25cm x 20cm)
- Ruthinium Fibra MAT - 10 sheets (10cm x 10cm)
- Ruthinium Fibra Base Resin - bottle (25ml)
- Ruthinium Fibra Catalyst - bottle (25ml)

### **RUTHINIUM FIBRA UHM – ULTRA HIGH MODULUS**

Instructions for use

#### Development of substructure

- Starting from a master mould with the abutments in position (duplication is not necessary), make a paper template not much larger than the substructure to be produced.
- Activate the resin at a ratio of 3:1, with 3 parts Ruthinium Fibra Base Resin and 1 part Ruthinium Fibra Catalyst.
- Completely impregnate the UHM fabric sheet with resin. Distribute the resin using the silicone spatula, ensuring that the fabric is completely impregnated. Remove any excess material.
- With a pair of scissors, trim the outer edges bordered by the tape. The fabric will not fray as it is impregnated.
- Place the paper template on the impregnated fabric and cut a rectangle large enough to contain the template.
- Cut the whole UHM fabric sheet as in the previous point to obtain more rectangles.
- Overlap all the rectangles obtained (a minimum of 13), changing by 10-15 degrees the direction of the fibres between each sheet.
- Compress the overlapped sheets using the roller provided (lamination).
- Place the paper template on the lamination sheets and trim them around the edge of the template. Place the obtained result on the master mould and extend with a pointed instrument near the abutments so as to fit the laminated sheets.
- Increase the volume of the substructure with the excess fibre from the laminated cut-off.
- Insert the master mould with the fibre placed inside a (food grade) vacuum bag and seal after vacuuming.
- Put in water at room temperature, bring to 80°C and keep the product at a temperature of 80°C for 2 hours.
- After the post-curing cycle, remove the product, allow it to cool and finish it using suitable cutters.
- Sandblast with 110/130 micron aluminium oxide. Brush the product with 96% ethanol. Important: do not vaporise.
- Apply a layer of light-curing primer and wait for it to be absorbed by the substructure.
- Apply a second layer of primer and blow to remove any excess. At this point it is possible to light-cure following the manufacturer's indications.
- Opacity with a light-curing matt, observing the manufacturer's instructions.
- At this point the denture can be constructed on top of the structure created.

#### Production of SPECIAL reinforcement

- Make a paper template not much larger than the substructure to be produced.
- Mould an elastic thermoplastic disk to the desired thickness. Use material that acrylic does not adhere to.
- Make some stops
- Limit the amount of fabric to be impregnated with paper tape, calculating 4 sheets to be laminated.
- Activate the resin at a ratio of 3:1, with 3 parts Ruthinium Fibra Base Resin and 1 part Ruthinium Fibra Catalyst.
- Impregnate the UHM fabric sheet part with resin. Distribute the resin using the silicone spatula, ensuring that the fabric is completely impregnated. Remove any excess material.
- With a pair of scissors, trim the outer edges bordered by the tape. The fabric will not fray as it is impregnated.
- Place the paper template on the impregnated UHM fabric and cut a rectangle large enough to contain the template.
- Cut the whole UHM fabric sheet as in the previous point to obtain 4 rectangles of fabric.
- Overlap all the rectangles obtained, changing by 10-15 degrees the direction of the fibres between each sheet.
- Compress the overlapped sheets using the roller provided (lamination).
- Place the paper template on the lamination sheets and trim them around the edge of the template.
- Fill the stops with opaque colour self-curing resin (powder), e.g. Acry Self OP1 colour, mixed with the excess resin used to impregnate the fabric.
- Place the laminated sheets obtained on the moulded thermoplastic disk.
- Insert the master mould with the fibre placed inside a (food grade) vacuum bag and seal after vacuuming.
- Put in water at room temperature, bring to 80°C and keep the product at a temperature of 80°C for 2 hours.
- After the post-curing cycle, remove the product, allow it to cool and finish it with suitable cutters.
- Sandblast with 110/130 micron aluminium oxide. Brush the product with 96% ethanol. Important: do not vaporise.
- Apply a layer of light-curing primer and wait for it to be absorbed by the substructure.
- Apply a second layer of primer and blow to remove any excess. At this point it is possible to light-cure following the manufacturer's indications.
- Opacity with a light-curing matt, observing the manufacturer's instructions.
- At this point the denture can be constructed on top of the structure created.

## RUTHINIUM FIBRA MAT

### Instructions for use

- Make a paper template not much larger than the substructure to be produced.
- Mould an elastic thermoplastic disk to the desired thickness. Use material that acrylic does not adhere to.
- Make some stops and insulate the mould with a silicone insulator at the stops.
- Place the previously made paper template on the fabric and make the reinforcement shape.
- Activate the resin at a ratio of 3:1, with 3 parts Ruthinium Fibra Base Resin and 1 part Ruthinium Fibra Catalyst.
- Distribute the resin, ensuring that the MAT fabric sheet is completely impregnated. Remove any excess material with absorbent paper.
- Fill the stops with opaque colour self-curing resin (powder), e.g. Acry Self OP1 colour, mixed with the excess resin used to impregnate the fabric.
- Place the laminated sheets obtained on the moulded thermoplastic disk.
- Insert the master mould with the fibre placed inside a (food grade) vacuum bag and seal after vacuuming.
- Put in water at room temperature, bring to 80°C and keep the product at a temperature of 80°C for 2 hours.
- After the post-curing cycle, remove the product, allow it to cool and finish it with suitable cutters.
- Sandblast with 110/130 micron aluminium oxide. Brush the product with 96% ethanol. Important: do not vaporise.
- Apply a layer of light-curing primer and wait for it to be absorbed by the substructure.
- Apply a second layer of primer and blow to remove any excess. At this point it is possible to light-cure following the manufacturer's indications.
- Opacify with a light-curing matt, observing the manufacturer's instructions.
- At this point the denture can be constructed on top of the structure created.

### CAUTION:

- Always observe a mixing ratio of 3:1 between Ruthinium Fibra Base Resin and Ruthinium Fibra Catalyst, regardless of the quantities used.
- Always perform the post-curing cycle according to the procedures previously described.
- Observe the following values in the design phase:
  - minimum thickness around the abutment: 1.5mm
  - minimum bar thickness: 5mm x 7mm
- For the substructures: maximum one tooth per cantilever.
- Do not leave any areas of fabric unimpregnated with epoxy resin.
- Do not deform the fibres or let them fray during processing.
- Do not cut the UHM fabric with scissors if not completely impregnated with the resin.
- Do not extend the UHM fibres too much so that the abutment enters the laminated fabric, in order to avoid lamination defects.
- Do not allow water to come into contact with the carbon fibre during post-curing.
- Post-curing must take place within an hour and a half after mixing the base and the catalyst.
- Always use disposable gloves when handling the epoxy resin. The disposable gloves must be powder-free to avoid contaminating the product.
- Do not vaporise the carbon fibre.
- Use gloves, a mask and a vacuum cleaner when finishing the fibre after post-curing.

### STORAGE

Ruthinium Fibra UHM and MAT: store in a clean and dry place.

Ruthinium Fibra Base Resin and Ruthinium Fibra Catalyst: store in a clean and dry place at a temperature between 15°C and 35°C.

The product is designed for professional use in the dental sector. Follow the instructions for use. The manufacturer shall accept no responsibility for damages arising from different or inadequate use. The user must personally check the suitability of the product for the intended use, particularly if these uses are not reported in the instructions for use. Do not use after expiration.

Warning: Any serious accident that may occur with the device must be reported by the user to the manufacturer and the competent authority of the Member State where the user and/or patient are established.

### Ruthinium Fibra Resina Base:

Causes skin irritation. Causes serious eye irritation. May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects. Avoid breathing vapours. Wear protective gloves / face protection. IF ON SKIN: Wash with plenty of water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention.  
Contains: reaction product: bisphenol-A-(epichlorhydrin) epoxy resin (number average molecular weight ≤700) Reaction products of hexane-1,6-diol with 2-(chloromethyl)oxirane (1:2)

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
Causes severe skin burns and eye damage. May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing vapours. Wash hands thoroughly after handling. Wear protective gloves / eye protection / face protection. IF ON SKIN: Wash with plenty of water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER/doctor if you feel unwell.  
Contains: 4,4'-Isopropylidenediphenol, oligomeric reaction products with 1-chloro-2,3-epoxypropane, reaction products with 3-aminoethyl-3,5,5-trimethylcyclohexylamine 3-aminomethyl-3,5,5-trimethylcyclohexylamine

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Ruthinium® - Dental Manufacturing S.p.A.