Ruthinium **DISC**

PMMA CAD-CAM DISC

Ruthinium **DISC**
Ruthinium **DISC** MULTILAYER
Ruthinium **DISC** MULTILAYER SC

Ruthinium **DISC** BiColor
Ruthinium **DISC** BiColor MULTILAYER

DISCO IN PMMA PER LA LAVORAZIONE CAD-CAM
PMMA DISC FOR CAD-CAM PROCESSING
DISCO PMMA PARA PROCESSAMENTO CAD-CAM
DISQUE PMMA POUR TRAITEMENT CAD-CAM
DISCO PMMA PARA LA ELABORACIÓN DE CAD-CAM
PMMA-SCHEIBE FÜR DIE CAD-CAM-BEARBEITUNG

(I) Istruzioni d'uso
(GB) Instructions for Use
(P) Instruções de uso
(F) Mode d'emploi
(E) Instrucciones de uso
(D) Gebrauchsanweisung
(NL) Gebruiksaanwijzing
(HU) Használati útmutató
(RO) Instrucțiuni de folosire
(TR) Kullanım için talimatlar
(GR) Οδηγίες χρήσης

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 Fabbricante

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(GB) PMMA DISC FOR CAD-CAM PROCESSING FOR THE MANUFACTURE OF TEMPORARY TEETH, ANATOMICAL BRIDGES AND/OR CROWNS, BASES FOR REMOVABLE PROSTHESES, DENTAL TEMPLATES AND ORTHODONTIC SPLINTS INTENDED TO COMPENSATE MASTICATORY HANDICAPS

INDICATIONS: Used for the fabrication of dental prostheses, including radiopaque ones. The disc is made of high molecular weight cross-linked resins that faithfully reproduce details while maintaining strength even for thin thicknesses. Can also be used for the fabrication of castable parts.

CONTRAINDICATIONS:

- Not indicated in cases of suspected allergy and/or proven sensitivity to the components of Ruthinium Disc.
- All uses not indicated.

AVAILABLE FORMATS:

WIEBER ATTACHMENT (Ø 98,5 mm) thicknesses: 14mm - 16mm - 18mm - 20mm - 25mm - 30mm / Ruthinium Disc Bicolor: thickness 30mm
 ZIRKOZAHN ATTACHMENT (Ø 95 mm) thicknesses: 16mm - 18mm - 20mm - 25mm - 30mm / Ruthinium Disc Bicolor: thickness 30mm
 AMMAN GIRRBACH ATTACHMENT (Ø 72 mm) thicknesses: 16mm - 20mm - 25mm - 30mm / Ruthinium Disc Bicolor: thickness 30mm

COLOURS AVAILABLE:

Ruthinium Disc: A1 - A2 - A3 - A3,5 - A4 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4 - PINK - NEUTRAL - RX - BL
 Ruthinium Disc Bicolor: A1-A2-A3-A3,5
 Ruthinium Disc Multilayer/ Ruthinium Disc Multilayer SC: A1 - A2 - A3 - A3,5 - A4 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4
 Ruthinium Disc Bicolor Multilayer: A1-A2-A3-A3,5 other thicknesses and colours available on request.

MAIN FEATURES:

- Medical device intended for dental CAD/CAM: the dental laboratory provides milling to obtain teeth, modelled structures such as crowns, bridges.
- High stability.
- Chromatic inalterability.
- Compliant with ISO 10477:2020

INSTRUCTIONS FOR USE:

1. The material must be processed in an environment with controlled conditions to avoid contamination of the material. In case of contamination by grease or oil, clean the material with alcohol.
2. Vacuum the chips during machining, avoid overheating the material to prevent deformation.
3. Refinish as usual, with burs, polishers, brushes and with pumice. Polish to a high gloss using brushes and polishing products (Acry Lux Crema, Acry Lux Liquido, Acry Lux Pasta), avoiding excessive overheating.

DESIGN SPECIFICATIONS:

- Minimum crown thickness: 3.5 mm occlusal and 0.6 mm circumferential;
- Offset: 3.0 mm;
- Distance to outer edge: 2.0 mm;
- Bridge connector radius: between 1.5 mm and 2.0 mm;
- Bridge connector cross-section: between 7.0 mm² 12.6 mm²

PERMITTED APPLICATIONS:

- Single elements: maximum 25 for molars, maximum 35 for incisors;
- Front and rear crowns;
- Front and rear bridges;
- Upper and lower arches;

WARNINGS:

- Before using the artefact in the oral cavity, proceed with cleaning and/or as defined by the doctor/dentist;
- The product is not autoclavable;
- Avoid cleaning with steam;
- All commercially available eugenol-free cements or adhesives can be used for cementation.

STORAGE CONDITIONS: We recommend keeping the product in a cool, dry place. PMMA Disc is not a hazardous product. Disposal must be in accordance with the laws of the competent authorities in force at the time of the activity. The waste can be treated as urban waste or recycled as normal plastic unless otherwise stated by the legislation in force at the time of the activity.

The product is designed for professional use in the dental sector. Observe the instructions in the instructions for use. The manufacturer accepts no liability for damage resulting from different or inappropriate use. All uses other than the intended use of the medical device are prohibited. Do not use after the expiry date.
 Warning: The user is required to report any accident and/or serious incident that may occur with the device both to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.



The physician/dentist should give instructions on the correct use of the medical device.

RECOMMENDED WORKING STRATEGY

TOOL TYPE							MILLING PARAMETERS				
Ø TOOL mm	Ø d mm	Nr. of Blades	R mm	l mm	i mm	L mm	RPM (S)	TOOL FEED XY (F)	TOOL FEED Z (F)	Ae (XY)	Ap (Z)
2	2	2	1	1.6	20	50	18500	1200	800	1	0.5
1	1	2	0.5	0.8	16	50	25000	1000	600	0.15	0.08
0,6	0,6	2	0.3	0,5	10	50	30000	500	500	0.05	0.05

ØD = Constrained by gripper measure of the milling machine

TOOL TYPE and MILLING PARAMETERS = recommended values and working strategy, consult the tool supplier.

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