# DOS

DIGITAL ANTONELLO CROCE OCCLUSAL SYSTEM

---

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### What is DACOS?

DACOS is a system designed for use in dental laboratories to facilitate the repositioning and mounting of 3D-printed dental models on analog articulators. DACOS employs an integrated magnetic system to ensure proper alignment of the printed models, maintaining their precise position.

### Intended Use:

DACOS is intended for dental laboratories utilizing CAD/CAM technologies for the design of removable or fixed dental prosthetics. The device is optimized for use with 3D-printed dental models created with SLA, LCD, or DLP technologies, enabling proper mounting of the models on the most commonly used analog articulators.

### Objective:

The DACOS system aims to streamline the articulation of 3D-printed models on analog articulators for inspection and prosthetic finalization purposes. It ensures precision and repeatability in aligning dental models through the use of specific attachments and magnetic plates (DACOS Plates), enabling accurate occlusal and functional verification of digitally produced models.

### Scope of Application:

DACOS is recommended for the remounting of removable, fixed dental prosthetics, and complex dentures, including:

- Complete or partial dentures
- Extensive bridges - Bite registrations and occlusal splints
- Implant-supported - dental structures

### Warnings and Limitations:

DACOS is exclusively designed for professional use in dental laboratories and is not intended for direct use on patients. Proper use of the DACOS system requires initial calibration of the articulator according to the manufacturer's instructions.

### Compatibility:

DACOS is compatible with the leading articulator systems through the self-installation mechanism of the DACOS plates.

## DACOS Hardware Installation Instructions

This guide provides detailed instructions for installing DACOS Plates in your articulator.

For software configuration and usage, refer to the video tutorials available on Blender For Denta: <https://www.blenderfordental.com/shop>.

### Required Materials:

- DACOS Plates ordered from Ruthinium.
- Pilot Calibrator (provided in the kit).
- Master Calibrator (downloadable as an STL file and 3D printable). 3D printer and model resin.
- Plaster + Silicone or wax for temporary fixation.
- GlueJet (included in the kit).
- Additional magnets or elastic bands (optional).
- Cyanoacrylate adhesive (for magnet fixation).

### 1 Printing the Master Calibrator

Before beginning the installation, print the Master Calibrator using your 3D printer and chosen model resin.

Select the appropriate size: Tight, Standard, or Loose, available within the DACOS software.

Ensure you use the same 3D printer and resin that will later be employed for printing your models. This ensures precise calibration of the entire DACOS system.

### 2 Installing the DACOS Plate with the Pilot Calibrator

If the articulator has a built-in occlusal plane (Camper Plane or Frankfurt Plane):

#### Positioning the Upper DACOS Plate:

- Place the Upper DACOS Plate in the articulator using the Pilot Calibrator to align it with the desired occlusal plane.
- Temporarily fix the DACOS Plate in position using a small amount of wax for stabilization.

#### Fixing with Plaster:

- Apply plaster to the Upper DACOS Plate, ensuring an even base for stable retention.
- Allow the plaster to fully set for approximately 30 minutes.
- Ensure the articulator's incisal pin is set to 0 and locked in place.

#### Removing the Pilot Calibrator:

Gently remove the Pilot Calibrator and replace it with the printed Master Calibrator.

#### Installing the Lower DACOS Plate:

Attach the Lower DACOS Plate to the Master Calibrator and complete the plastering process, ensuring perfect alignment between the plates.

### 3 Installing the DACOS Plate Using Elastic Bands and the Bonwill Triangle

If the articulator does not have a built-in occlusal plane:

#### Preparing the Master Calibrator:

- Print the Master Calibrator (Tight, Standard, or Loose)
- Insert magnets into the calibrator's attachments using the GlueJet tool.
- Optionally, add additional magnets to the Master Calibrator for enhanced stability, or secure it with elastic bands or cyanoacrylate adhesive.

#### Aligning the Occlusal Plane:

Loop an elastic band around the articulator, connecting the two posterior notches and the anterior notch, forming the Bonwill Triangle. Use the Master Calibrator as a reference for the occlusal plane.

#### Temporarily Fixing the Lower DACOS Plate:

Place the Lower DACOS Plate and temporarily secure it with silicone, ensuring alignment with the Master Calibrator and the elastic plane.

#### Installing the Upper DACOS Plate:

- Ensure the articulator's incisal pin is set to 0 and locked in place.
- Proceed with plastering the Upper DACOS Plate.
- Once the plaster sets, remove the silicone and finalize the installation by plastering the Lower DACOS Plate.
- Allow approximately 30 minutes for the plaster to fully set.

### 4 Attaching Magnets to Models Using GlueJet

Once your models are generated using the DACOS module, follow these steps for finalizing the hardware installation:

#### Preparing the Models:

- Print the models using your desired 3D printer and resin.
- Clean the models thoroughly before proceeding.

#### Placing the Magnets:

- Attract the Magnet: Use the GlueJet tool to pick up a magnet from the kit. The GlueJet tool automatically ensures the correct polarity for alignment with the DACOS Plates.
- Apply a small amount of cyanoacrylate adhesive inside the attachment slot on the model.

#### Inserting the Magnet:

- Insert the GlueJet tool with the magnet into the attachment slot.
- Remove the GlueJet Plunger, allowing the magnet to drop into the correct position.
- Repeat this process for all four attachment points on the model. Allow the adhesive to cure for approximately 10 minutes.

### 5 Conclusion

The models are now ready to be placed in the DACOS Plate installed in your articulator. Once positioned, you can proceed with occlusal checks directly in the articulator.

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