Ruthinium **DISC**

PMMA CAD-CAM DISC

Ruthinium **DISC** Ruthinium **DISC** MULTILAYER Ruthinium **DISC** MULTILAYER SC

DISCO IN PMMA PER LA LAVORAZIONE CAD-CAM
PMMA DISC FOR CAD-CAM PROCESSING
DISCO PMMA PARA PROCESSAMENTO CAD-CAM
DISQUE PMMA POUR TRAITEMENT CAD-CAM
DISCO PMMA PARA LA ELABORACIÓN DE CAD-CAM
PMMA-SCHEIBE FÜR DIE CAD-CAM-BEARBEITUNG

(I) Istruzioni d'uso
(GB) Instructions for Use
(P) Instruções de uso
(F) Mode d'emploi
(E) Instrucciones de uso
(D) Gebrauchsanweisung
(NL) Gebruiksaanwijzing
(HU) Használati útmutató
(RO) Instrucțiuni de folosire
(TR) Kullanım için talimatlar
(GR) Οδηγίες χρήσης

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

Via Cà Mignola Nuova 1699 - 45021 Badia Polesine (RO) ITALY

Tel. +39 0425 51628 - www.ruthinium.it - info@ruthinium.it

INTENDED USE: PMMA DISC FOR CAD-CAM PROCESSING. INTENDED FOR THE MANUFACTURE OF DEVICES DESIGNED TO COMPENSATE FOR MASTICATORY DYSFUNCTIONS SUCH AS: TEETH FOR REMOVABLE AND NON-REMOVABLE DENTURES, NON-REMOVABLE ANATOMICAL BRIDGES AND/OR CROWNS, BASES FOR REMOVABLE AND NON-REMOVABLE DENTURES, ORTHODONTIC BITE SPLINTS. THE RADIOPAQUE PMMA DISC IS INTENDED TO PROVIDE A RADIOGRAPHICALLY VISIBLE SOLUTION TO FACILITATE DIAGNOSIS AND POST-TREATMENT CLINICAL MONITORING.

INSTRUCTIONS FOR USE:

The indication for the use of the PMMA CAD/CAM disc includes cases of missing, damaged, or compromised teeth where intervention is necessary to maintain functionality and aesthetics until the final treatment can be completed. The PMMA Disc is indicated for the fabrication of radiopaque dentures and radiological investigations.

EXPECTED CLINICAL BENEFITS:

Restoration of masticatory, phonetic, and aesthetic function through dentures
Orthodontic correction of masticatory dysfunctions and precision in radiological investigations

TARGET POPULATION:

PMMA discs are intended for all age groups and genders. Patients are individuals with masticatory dysfunctions requiring interventions to restore chewing through the clinical application of dentures, partial or complete anatomical arches, bridges and crowns, bases for removable dentures, surgical guides, and orthodontic bite splints.

PROFILE OF PROFESSIONAL USERS:

The recipients of the PMMA Disc are dental professionals and technicians with specific skills in digital dentistry, CAD-CAM technology, and the fabrication of dentures, teeth, anatomical bridges, crowns, bases for removable dentures, dental masks, and orthodontic bite splints.

CONTRAINDICATIONS:

Not indicated in case of suspected allergy and/or proven sensitivity to components of the Ruthinium Disc. All uses not specified in the instructions for use.

AVAILABLE FORMATS:

WIEBER ATTACHMENT (Ø 98.5 mm) thicknesses: 14mm - 16mm - 18mm - 20mm - 25mm - 30mm

ZIRKOZAHN ATTACHMENT (Ø 95 mm) thicknesses: 16mm - 18mm - 20mm - 25mm - 30mm

AMMAN GIRRBACH ATTACHMENT (Ø 72 mm) thicknesses: 16mm - 20mm - 25mm - 30mm

AVAILABLE SHADES:

Ruthinium Disc: A1 - A2 - A3 - A3,5 - A4 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4 - PINK - NEUTRAL - RX - BL - BLUE - ORNG - YLLW - ROSE

Ruthinium Disc Multilayer/ Ruthinium Disc Multilayer SC: A1 - A2 - A3 - A3,5 - A4 - B1 - B2 - B3 - B4 - C1 - C2 - C3 - C4 - D2 - D3 - D4

MAIN FEATURES:

- Medical device intended for dental CAD/CAM: the dental laboratory provides milling to obtain teeth, modelled structures such as crowns, bridges.
- High stability and chromatic inalterability. Compliant with ISO 10477:2020

INSTRUCTIONS FOR USE:

1. The material must be processed in an environment with controlled conditions to avoid contamination of the material. In case of contamination by grease or oil, clean the material with alcohol.
 2. Vacuum the chips during machining, avoid overheating the material to prevent deformation.
 3. Refinish as usual, with burs, polishers, brushes and with pumice. Polish to a high gloss using brushes and polishing products (Acry Lux Crema, Acry Lux Liquido, Acry Lux Pasta), avoiding excessive overheating.
- The device might be used for fixing with: •Temporary cements/fixing materials, eugenol-free •Temporary cementation •Permanent adhesive
- The fixing of Ruthinium Disc must always be carried out according to the instructions for use of the respective manufacturer.

DESIGN SPECIFICATIONS:

- Minimum crown thickness: 3.5 mm occlusal and 0.6 mm circumferential;
- Offset CAM/Machine route: 3.0 mm;
- Minimum distance to outer edge: 2.0 mm;
- Bridge connector radius: between 1.5 mm and 2.0 mm;
- Radius of connectors for surgical templates and plates 3.0 mm and 4.0 mm;
- Bridge connector cross-section: between 7.0 mm² 12.6 mm²
- Cross section connectors for surgical guides and basal plates 14 mm² and 25.2 mm²

CHEMICAL COMPOSITION:

Polymethyl methacrylate	> 98.5%
Methyl methacrylate	< 1%
Color pigments	< 0.5%

PERMITTED APPLICATIONS:

- Single elements: maximum 25 for molars, maximum 35 for incisors;
- Front and rear crowns; Front and rear bridges; Upper and lower arches;

WARNINGS:

- Before using the artefact in the oral cavity, proceed with cleaning and/or as defined by the doctor/dentist;
- The product is not autoclavable; Avoid cleaning with steam; All commercially available eugenol-free cements or adhesives can be used for cementation.

STORAGE CONDITIONS: We recommend keeping the product in a cool, dry place. PMMA Disc is not a hazardous product. Disposal must be in accordance with the laws of the competent authorities in force at the time of the activity. The waste can be treated as urban waste or recycled as normal plastic unless otherwise stated by the legislation in force at the time of the activity.

The product is designed for professional use in the dental sector. Observe the instructions in the instructions for use. The manufacturer accepts no liability for damage resulting from different or inappropriate use. All uses other than the intended use of the medical device are prohibited. Do not use after the expiry date.

Warning: The user is required to report any accident and/or serious incident that may occur with the device both to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established. The physician/dentist should give instructions on the correct use of the medical device.



For further information on the safety and clinical performance of this device, please consult the Summary of Safety and Clinical Performance (SSCP) available in the European EUDAMED database:<https://ec.europa.eu/tools/eudamed>

RECOMMENDED WORKING STRATEGY

TOOL TYPE							MILLING PARAMETERS				
Ø TOOL mm	Ø d mm	Nr. of Blades	R mm	l mm	i mm	L mm	RPM (S)	TOOL FEED XY (F)	TOOL FEED Z (F)	Ae (XY)	Ap (Z)
2	2	2	1	1.6	20	50	18500	1200	800	1	0.5
1	1	2	0.5	0.8	16	50	25000	1000	600	0.15	0.08
0,6	0,6	2	0.3	0,5	10	50	30000	500	500	0.05	0.05

ØD = Constrained by clamps measurement of the milling machine

TOOL TYPE and MILLING PARAMETERS = recommended values and working strategy, consult the tool supplier.

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Ruthinium® - Dental Manufacturing S.p.A.